

RNP Nucleofection for cell lines using Lonza 4D Nucleofector

Protocol:

1. For each 20ul nucleofection aliquot 200k cells into 15 mL conical tubes.
2. Prepare RNP mix. Cas9 buffer is kept in the TC hood and must be kept sterile.
 - a. Bring 100 pmol of Cas9 (15ug) to a final volume of 5 μ L using Cas9 buffer (20 mM HEPES-KOH pH 7.5, 150 mM KCl, 10% glycerol, 1 mM TCEP). Add TCEP fresh each time.
 - b. Bring 120 pmol sgRNA (~4ug) to a final volume of 5 μ L using Cas9 buffer.
 - c. Add Cas9 to sgRNA very, very slowly while swirling pipette tip, should take 30 s to 1 minute.
 - d. Allow RNP complex to form at room temp for 10-20 minutes. During this step move on to preparing cells below.
3. Prepare Cells
 - a. Spin cells at 100 x g (or below) for 10 minutes to pellet cells softly. While the cells are spinning, prepare plate and cuvette.
 - b. Prepare a 12-well-plate with 1mL media per well, and pre-warm in the incubator.
4. Nucleofection
 - a. Prepare and label wells on 20uL nucleofection strips. Configure Lonza 4d using recommended cell-type program.
 - b. Pipette off media from cells, gently but completely, using a P200. The pellet is very soft so be careful.
 - c. Resuspend cells in 20 μ L of nucleofector solution using a P200.
 - d. Add the entire 10 μ L RNP mix to the 20 μ L resuspension and mix.
 - e. Optional: Add 1uL of 100uM donor DNA (100 pmoles) and mix well.
 - f. Add nucleofection mixes to the multiwell cuvette, and cap.
 - a. Pay attention to the orientation of the cap and cuvette in the nucleofector, which is noted in the manufacturer's instructions.
 - b. Add carefully to one short side of the well, at an angle.
 - c. Do not produce any bubbles. The solution does not need to be completely filling the well as long as there are no bubbles.
 - g. Insert cuvette into nucleofector and run nucleofection program.
 - h. Allow cells to sit in nucleofection strips for 10 minutes post-nucleofection. This is supposed to increase efficiency.
 - i. Add 80uL of pre-warmed media to each well.
 - j. Pipette mixture out with a P200 into your pre-warmed 12-well plate. This should get the vast majority of cells, but if you wish, you may wash out the rest with media from the same well, chemistry-style.
 - k. Allow cells 24 hours to settle and recover before moving on to downstream analysis.