

HSC Room 60:

NIS elements version 6.02 (one floating license)

- Nikon Image Analysis software suite
- View, extract, stitch, pre-process Nikon “nd2” files
- NIS denoise.ai for deep learning denoising
- General Analysis 3 (GA3) module for building image analysis workflows (processing, segmentation, measurement, data management, etc.)

HCI Room 1470:

NIS elements version 6.02 (one floating license)

- Nikon Image Analysis software suite
- View, extract, stitch, pre-process Nikon “nd2” files
- NIS denoise.ai for deep learning denoising
- General Analysis 3 (GA3) module for building image analysis workflows (processing, segmentation, measurement, data management, etc.)

ZEN Desk version 3.8 (local license)

- Zeiss Image Analysis software suite
- With 2D Toolkit and Bio Apps modules for creating 2D image analysis pipelines
- Deconvolution module (basic)
- Macro processing (e.g. for manual spectral unmixing / background subtraction)